

DiNovo

Pharmacy & Packaging Provisions



Distributing Innovation and Compliance



www.agapilabsolutions.com

Contact Us:

Phone 1 (267) 397-3345

E-mail customerservice@agapilabsolutions.com

DiNovo

Pharmacy & Packaging Provisions

The NEW nationwide distributor of USP <797> and <71> compliance products, packaging solutions, and medical supplies designed to bring innovative ideas at a low cost to the healthcare market place.

DiNovo's goal, unlike other medical distributors, is not to move boxes but to become a valued partner and resource to help aid in pharmacy compliance. We are proud to represent some of the top manufacturers in the medical industry such as Medline, ITW Texwipe, Hardy Diagnostics, and many more.

Our core product offering is as followed but not limited to:

- Adhesive Mats & Frames
- Sterile Coveralls
- Protective Apparel
- Disposal Microfiber Mop
- Reusable PVA Roller Mop
- Employee Validation Kits
- Environmental Monitoring
- Sterility Testing
- Sterile Empty Vials
- Sterile Alcohol
- Surface Disinfectants
- Dry/Wet Cleanroom Wipers
- Sterile Tamper Evident Seals
- Tamper Evident Caps
- Custom Polyethylene Shippers
- Cold Chain/Packaging Products

Adhesive Mats & Frames



CleanStep™

Adhesive Mats & Frames

Features

Mats

- Adhesive-coated, removable sheets made from high grade polyethylene
- Consistent adhesive coating from sheet to sheet
- Available in 5 sizes, 3 colors and 30 or 60-layer
- Numbered tabs take the guesswork out of inventory and replacement

Frames

- Portable polystyrene
- Non-slip design
- Lightweight
- Available in 3 sizes and 2 colors
- Protects floor surfaces

Description

CleanStep™ is the preferred choice in adhesive mats and frames for cleanrooms and critical areas. Constructed from the highest grade polyethylene (mats) and polystyrene (frames), this contamination control system grabs and retains particulates from foot and wheel-borne contaminants.

Adhesive Mats

Each removable, adhesive-coated polyethylene sheet traps particles and prevents transfer. CleanStep™ mats feature an antimicrobial agent to prevent bacterial growth and numbered tabs to simplify inventory and replacement.

Frames

These light-weight, non-skid polystyrene frames adhere to the floor without adhesive and stay in place utilizing a high-grade non-skid backing. Easy to move for applications needing portability. Frames fit 18" x 36", 25" x 45" mat sizes. Available in black or white color choices.

Applications

- Cleanrooms
- Critical environment entry points
- Gowning rooms
- Ante rooms
- Cleanroom construction projects
- Laboratories

Product	Item	Description	Pack
Adhesive Mats	AMA183682W	CleanStep, Tacky Mat, White, 18x36	60 layers/mat, 4 mats/case
	AMA183682B	CleanStep, Tacky Mat, Blue, 18x36	60 layers/mat, 4 mats/case
	AMA254582W	CleanStep, Tacky Mat, White, 25x45	60 layers/mat, 4 mats/case
	AMA254582B	CleanStep, Tacky Mat, Blue, 25x45	60 layers/mat, 4 mats/case
Frames	FB1937	Frame, Black 19x37	1 Frame
	FW1937	Frame, White 19x37	1 Frame
	FB2646	Frame, Black 26x46	1 Frame
	FW2646	Frame, White 26x46	1 Frame



www.agapilabsolutions.com

Contact Us:

Phone 1 (267) 397-3345

E-mail customerservice@agapilabsolutions.com

Sterile Protective Coveralls

ProGard offers a comfortable fabric with a superior liquid barrier, low linting and strong tear resistance. All ProGard products are manufactured in an ISO-certified facility. They are sterilized to a 10⁻⁶ Sterility Assurance Level and individually packaged and sealed for next-level protection.

ProGard Products are made with serged seams and tunnelized elastic to strengthen the seam, provide secure closure and eliminate raw edges that may generate particulates.



Product Features

- Attached hood and boots
- Serged seams
- Anti-skid resistant boots
- Elastic wrist
- Zipper with storm-fly
- Closure at chin

Key Benefits

- Sterility Assurance 10⁻⁶
- Extreme low linting
- Recyclable
- More Breathable
- Liquid barrier
- Economical

Product	Item	Description	Pack
Sterile Protective Coveralls	DYNDA1229	Bunny Suit, Sterile, Individual Wrapped, Over Pack, Medium	25/cs
	DYNDA1230	Bunny Suit, Sterile, Individual Wrapped, Over Pack, Large	25/cs
	DYNDA1228	Bunny Suit, Sterile, Individual Wrapped, Over Pack, X-Large	25/cs
	DYNDA1227	Bunny Suit, Sterile, Individual Wrapped, Over Pack, XX- Large	25/cs
	DYNDA1226	Bunny Suit, Sterile, Individual Wrapped, Over Pack, XXX-Large	25/cs




www.agapilabsolutions.com

Contact Us:

Phone 1 (267) 397-3345

E-mail customerservice@agapilabsolutions.com



Disposal Protective Apparel

We know you rely on disposable apparel to keep your staff protected and comfortable. For this reason, DiNovo provides products with the finest materials available for protection, offering the greatest level of comfort.

- Bouffant Caps ● Shoe Covers
- Gowns ● Face Masks ● Gloves
- Sleeve Covers ● Beard Covers



Product	Item	Description	Pack
Bouffant Caps	CRI1002	Bouffant Cap, 21", Blue	500/cs
	CRI1004	Bouffant Cap, 24", Blue	500/cs
Shoe Covers	CRI2000	Shoe Cover, Breathable, Spunbond, Fluid Repellent, Regular	300/cs
	CRI2001	Shoe Cover, Breathable, Spunbond, Fluid Repellent, X-Large	200/cs
Gowns	CRI4000	Gown, Spunbond, Fluid Repellent, Waist-Neck Tie, Yellow, Regular	50/cs
	CRI4001	Gown, Spunbond, Fluid Repellent, Waist-Neck Tie, Yellow, X-Large	50/cs
Face Masks	NON27122	Face Mask, Yellow, Ear loop	300/cs
	NON27376	Face Mask, Blue, Ties	300/cs
Gloves	Sterile	We carry a full range of Sterile Nitrile, Vinyl, and Latex Gloves	
	Non-Sterile	We carry a full range of NonSterile Nitrile, Vinyl and Latex Gloves	
Sleeves	DYNJP2000	Sleeve, Sterile, 23" Length, Individual	60/cs
Beard	NONSH400	Beard Cover, Spunbond, Elastic, White	1000/cs

DiNovo is a full service Medline Distributor and carries a large array of products. Contact DiNovo for all requests.



www.agapilabsolutions.com

Contact Us:

Phone 1 (267) 397-3345

E-mail customerservice@agapilabsolutions.com

Microfiber Mopping System

Maintaining floors, walls, and ceilings in critical environments have never been easier.

The NEW DiNovo microfiber mopping system is becoming the market standard for cleanroom mopping systems. Proper cleanroom mopping is extremely important while trying to maintain a critical environment.

The microfiber disposable mop system provides easily changeable mop pads. This arrangement permits users to rapidly change mop pads with different disinfectants, adhering to strict wiping protocol. The Microfiber disposable mop provides the cleanest substrate available for cleaning walls and ceilings. The disposable microfiber mop pads conforms to textured floors and other surfaces for the most efficient cleaning.

Features

- Microfiber mop pads are low-lint and extremely absorbent
- Can be used dry to attract dust and particles or wet for mopping
- The mop pads are disposable
- 60" telescopic mop pole, with cushion handle
- 16" swivel mop frame utilizes a locking device
- Flag System on handle and mop frame



Benefits

- Removes and disinfects bioburden from cleanroom surfaces
- Can be used with various disinfectants
- Protects against cross contamination
- Maintains employee safety by not over handling dirty mop pads
- Makes for a more comfortable cleaning experience for employees of all heights
- The locking device makes cleaning the walls and ceiling easier
- Can be used to help identify what mop is used in specific areas and/or with certain disinfectants

Product	Item	Description	Pack
Microfiber Mop System	DMH60	Mop Handle, Telescopic Up to 60", Flag System	1 Handle
	DMHF60	Mop Handle, Flow Mop , Bucketless System	1 Handle
	DMF16	Mop Frame, 16", Locking Device, Flag System	1 Mop Frame
	DMMP18	Microfiber Mop Pads, Disposable, 18", White	100 Pads/Case
Disinfectants	ASO646108Z	Vesphene, Germicidal, High PH Cleaner, Gallon	1 Bottle
	ASO646608Z	LpH, Germicidal, Low PH Cleaner, Gallon	1 Bottle
	TX6466H	TbC Tablets (Bleach Alternative), 1000 ppm chlorine	15/Bag

PVA Sponge Mopping System



The Super® PVA Mop features 2 sets of rugged, double roller wringers. Made for continuous, daily use by professionals and consumers. Not affected by most chemicals. Soft, absorbent, horizontal ribbed PVA sponge lifts dirt and debris. Static - free, lint - free and drip-free, the mop is lightweight and durable.

- The Super® mop heads are only made with 100% pure PVA sponge
- Absorbs chemical or liquid spills fast
- Applies disinfectant or solutions evenly on walls, ceilings without dripping
- Attracts lint, dirt, and other particulates like a magnet
- Patented, durable double wringers completely cover sponge when wringing to provide a clean, drip-free mop
- Adjustable aluminum shaft minimizes hand and arm exertion, extends from 39" to 52"; Optional Super Extension Handle extends from 50" to 73"
- Mop body and wringer handle are made of tough, virgin plastic. Corrosion-resistant screws secure sponge head to mop holder

Durable Double Wringers

The 4 rollers completely cover the sponge when wringing, to provide a clean, drip-free mop.




5.2 GALLON BUCKET

High impact polypropylene plastic for durability with large carry handle. Includes removable dirt screen.



Soft, Ribbed PVA Sponge

Patented grooves lift dirt and debris. Pure PVA sponge does not promote the growth of bacteria.



Product	Item	Description	Pack
PVA Mopping System	335-0022	16" PVA Sponge Mop Set (Mop Handle w/ 1 Refill)	1/ea
	335-0012	16" PVA Sponge Mop Replacement Head Refill	1/ea
	987-0010	5.2 Gallon Bucket with Removable Dirt Screen	1/ea
	335-0018	Mop Extension, Add up to 20"	1/ea
Disinfectants	ASO646108Z	Vesphene, Germicidal, High PH Cleaner, Gallon	1 Bottle
	ASO646608Z	LpH, Germicidal, Low PH Cleaner, Gallon	1 Bottle
	TX6466H	TbC Tablets (Bleach Alternative), 1000 ppm chlorine	15/Bag



www.agapilabsolutions.com

Contact Us:

Phone 1 (267) 397-3345

E-mail customerservice@agapilabsolutions.com

Employee Validation Kits

Microbial Contamination Testing for CSPs

DiNovo offers all the products you need to easily assess the risk of microbial contamination of your CSPs (Compounded Sterile Preparations) according to the USP Chapter <797> requirements for media-fill challenge testing.

- One kit supplies what is needed to test the proficiency of one technician or pharmacist
- Kits are available for low, medium, and high risk procedures
- Components are in a self contained ready-to-use kit
- Results log sheet is included
- Detailed instructions and diagrams are included
- Expert technical support is available from our staff of microbiologists
- Kits are designed to help you follow USP <797> guidelines

Low-Risk



Cat. no. HVL1

Medium-Risk



Cat. no. HVM1

Medium-Risk Basic



Cat. no. HVM2

High-Risk



Cat. no. HVH1

Compliance was never so easy!

Product	Item	Description	Pack
Low-Risk	HVL1	Low-Risk Level, CSP Challenge Kit	1/ea
Medium-Risk	HVM1	Medium-Risk Level, Comprehensive, CSP Challenge Kit	1/ea
	HVM2	Medium-Risk Level, Vial System, CSP Challenge Kit	1/ea
High-Risk	HVH1	High-Risk Level, CSP Challenge Kit	1/ea



www.agopilabsolutions.com


Contact Us:

Phone 1 (267) 397-3345

E-mail customerservice@agopilabsolutions.com

Environmental Monitoring

DiNovo offers a fast and convenient way to monitor the microbial bioburden of the work surface of your laminar air flow workbench. Contact agar plates are simply pressed onto the surface or glove fingertips, incubated, and colonies are counted.



Step 2

Surface Sampling




Gently press plate to surface without twisting or sliding

Finger Tip Sampling



Gently press glove tips to plate surface

Step 3



Replace lid and incubate at 33°C to 37°C. After 48 hours, count the number of colonies

Step 1



Remove lid

Product	Item	Description	Pack
Contact Plate	P34	Plate, Tryptic Soy Agar with Lecithin and Tween 80, 15x60mm	10/pk
Air Sampling	G60	Plate, Tryptic Soy Agar, 15x100mm	10/pk
	W28	Plate, Malt Extract Agar (selective for yeast and mold), 15x100mm	10/pk

DiNovo is an exclusive Hardy Distributor and sell an array of agar products. Contact us for all special requests.



www.agapilabsolutions.com

Contact Us:

Phone 1 (267) 397-3345

E-mail customerservice@agapilabsolutions.com

Sterility Testing

Direct Inoculation



Drug solutions can be tested for sterility by direct inoculation into Tryptic Soy Broth and Fluid Thioglycollate according to procedures found in <USP 71>. By incubating for 14 days, a lack of fungal, aerobic and anaerobic bacterial contamination can be proven.

Please consult <USP 71> (Table II) to verify proper inoculum size and media volume.

Below are some types of media that can be used for direct inoculation.

Tryptic Soy Broth, USP

Recommended for use as a general purpose medium for the isolation and cultivation of a wide variety of bacteria and fungi. 100ml glass bottle with needle port septum



Cat. No. U46

Fluid Thioglycollate Broth Medium with Indicator, USP

For the cultivation of anaerobic microorganisms. 100ml glass bottle with needle port septum



Cat. No. U84

Tryptic Soy Broth, USP


Recommended for use as a general purpose medium for the isolation and cultivation of a wide variety of bacteria and fungi. 20ml glass bottle with needle port septum



Cat. No. U38

Fluid Thioglycollate Broth Medium with Indicator, USP

For the cultivation of anaerobic microorganisms. 20ml glass bottle with needle port septum



Cat. No. U68

Product	Item	Description	Pack
TSB	U46	Tryptic Soy Broth, 100ml Glass Bottle	20/bx
	U38	Tryptic Soy Broth, 20ml Glass Bottle	50/bx
FTM	U84	Fluid Thioglycollate, Resazurin Indicator, Clear Formula, 100ml Glass Btl	20/bx
	U68	Fluid Thioglycollate, Resazurin Indicator, Clear Formula, 20ml Glass Btl	50/bx

DiNovo is an exclusive Hardy Distributor and can create custom containers and fill volumes upon request.



www.agopilabsolutions.com

Contact Us:

Phone 1 (267) 397-3345

E-mail customerservice@agopilabsolutions.com

Sterile Empty Vials



Empty Injection Vials are sterile and nonpyrogenic. These vials are made with USP Type 1 glass used for mixing different medications or solutions for injection.

Features

- Aluminum sealed top
- Butyl compound (synthetic rubber) stopper
- Amber Vials available in select sizes.
- Available with Color Rings: Blue, Yellow, Green, Orange, Purple, Red, Pink.



Product	Item	Vial Size	Stopper Size	Available In Amber	Available w/ Color Rings	Pack
Sterile Empty Vials	SV111	1ml	11mm			25/pk
	SV213	2ml	13mm	Yes	Yes	25/pk
	SV313	3ml	13mm			25/pk
	SV513	5ml	13mm		Yes	25/pk
	SV520	5ml	20mm	Yes	Yes	25/pk
	SV1013	10ml	13mm		Yes	25/pk
	SV1020	10ml	20mm	Yes	Yes	25/pk
	SV20	20ml	20mm			25/pk
	SV30	30ml	20mm	Yes		25/pk
	SV50	50ml	20mm	Yes		25/pk
SV100	100ml	20mm	Yes		25/pk	



www.agapilabsolutions.com

Contact Us:

Phone 1 (267) 397-3345

E-mail customerservice@agapilabsolutions.com

Sterile Isopropanol Alcohol

Sterile 70% isopropanol contains 70% by volume USP-grade isopropanol (isopropyl alcohol) and 30% USP-purified water. This solution is submicron filtered, filled into cleaned containers, double-bagged and gamma-irradiated. Each bottle is lot coded with an expiration date for easy recordkeeping.



Applications

- Surface cleaning
- Use on gloved hands in sterile suites
- Wipe down for pass through to controlled environment
- Ideally suited for use with sterile cleanroom wipers

Benefits

- Removes many surface contaminants
- Evaporates and leaves extremely low residues
- Prefiltered and sterilized for pass through and use in critical environments
- Eliminates mixing, filtration, sterilization, processing and QC expense of in-house production of sterile alcohol/water blends
- Ready to use sprayer eliminates handling and assembly of sprayhead and maintains sterility
- Meets FDA guideline for medical device bacterial endotoxin limit of not more than 20 EU/device

Features

- Made with USP-grade Isopropanol
- Filtered through 0.2 micron filter
- Gamma-irradiated to a sterility level of 10^{-6} with independent QC audits to ensure sterility
- Double-bagged in solvent-safe bags
- Non-atomizing adjustable spray
- Fully assembled, ready to use applicator / container
- Each bottle includes a set of 5 color indicator dots for easy identification
- Lot traceable
- Each lot tested for endotoxins
- Lot-specific information on each shipment simplifies recordkeeping:
 - »Certificate of Processing confirming radiation dosage.
 - »Certificates of Compliance and Analysis
 - »Expiration date marked on each package

Product	Item	Description	Pack
Sterile Alcohol	TX3290	Sterile 70% Isopropanol Alcohol, 1 gallon Polybottle	4 polybottles/case
	TX3273	Sterile 70% Isopropanol Alcohol, 32 fluid ounces Trigger Spray	12 polybottles/case
	TX3270	Sterile 70% Isopropanol Alcohol, 16 fluid ounces Trigger Spray	12 polybottles/case
	TX8270	Sterile 70% Isopropanol Alcohol, 8 fluid ounces Trigger Spray	12 polybottles/case



www.agopilabsolutions.com

Contact Us:

Phone 1 (267) 397-3345

E-mail customerservice@agopilabsolutions.com

Applications

Hard, non-porous surfaces

- Carts ● Chrome ● Non-Critical Equipment
- Glass ● Plastic ● Plexiglas® ● Polycarbonate
- Stainless Steel ● Racks/Shelves ● Windows

Surface Disinfectants

Quaternary Ammonium Disinfectant

TexQ® Disinfectant is the last generation of quaternary ammonium compounds (QACs), EPA registered one-step cleaner/disinfectant. It is effective against a broad spectrum of bacteria, viruses and fungi, and inhibits the growth of mold and mildew and their odors when used as directed. The bottles are gamma-irradiated so it is ready to be used in your cleanroom upon opening with functional use label on bottle for easy documentation and usage control. TexQ is free of dyes and fragrance so to protect employees of respiratory irritation and is pH neutral making it safe for operation and not aggressive on surfaces.



Product	Item	Description	Pack
Surface Disinfectant	TX650	Quaternary Ammonium, Ready-to-Use, 22oz. Trigger	1/ea
	TX651	Quaternary Ammonium, Concentrate, Gallon Bottle	1/ea

Bleach Germicidal Cleaner

Clorox Healthcare® Bleach Germicidal Cleaners are unique, stabilized sodium hypochlorite and detergent solutions that kill pathogens fast — with a one-minute contact time for a broad range of microorganisms, including MDROs such as *Acinetobacter baumannii*, *Enterobacter aerogenes*, VRE and MRSA.

Clorox Healthcare® Bleach Germicidal Cleaners are premixed, ready to use and are similar to the 1:10 bleach solution strength that is recommended by CDC guidelines for effective disinfecting. Clorox Healthcare® Bleach Germicidal Cleaners contain an anticorrosive ingredient that makes them less damaging to surfaces or equipment than conventionally prepared bleach solutions.



Product	Item	Description	Pack
Surface Disinfectant	CL-68970	Bleach Germicidal Cleaner, 32oz. Trigger Spray	1/ea
	CL-35309	Bleach Germicidal Cleaner, Wipes	420/cs



www.agapilabsolutions.com

Contact Us:

Phone 1 (267) 397-3345

E-mail customerservice@agapilabsolutions.com

Dry Cleanroom Wipers

Description

TexVantage® nonwoven wipers contain a blend of 65% cellulose (for sorption capacity) and 35% polyester (for strength and low particle release) for use in spill control and general-purpose cleaning.

Application

- ISO Class 5-8
- Aqueous spill cleanup
- General-purpose wiping
- Applying common cleaning and disinfecting solutions
- Ideal for wipe downs in sterile gowning rooms

Features

- Nonwoven, hydroentangled construction with bidirectional strength
- Solvent-safe cleanroom packaging
- Chemically resistant to withstand exposure to solutions such as bleach, phenols and quats
- Statistical process control of wiper quality
- Autoclavable
- Double-Bagged



TexVantage™

Polyester/Cellulose Wipers

Product	Item	Description	Pack
65% Cellulose	TX8939	TexVantage™, 9x9	300 wipes/bag, 20 bags/cs
35% Polyester	TX8932	TexVantage™, 12x12	150 wipes/bag, 20 bags/cs

DiNovo is an exclusive ITW Texwipe Distributor and sell an array of dry wipers. Contact us for all special requests.

Description

TechniCloth® wipers are fabricated from a hydroentangled, nonwoven blend of 55% cellulose and 45% polyester. This fabric blend combines the highly absorbent properties of a natural fiber with the cleanliness and strength of a synthetic wiper.

Application

- Designed for general-purpose wiping where economy and cleanliness are of critical importance
- Ideal for applications where high absorbency and low extractables are prime considerations
- Excellent for picking up aqueous spills
- Suitable for cleaning precision components and lab apparatus

Features

- Excellent wet and dry strength due to the polyester
- Highly absorbent due to the cellulose
- Nonwoven, hydroentangled construction with excellent bidirectional strength
- Low extractable levels
- Autoclavable
- Double-Bagged



TechniCloth®

Wipers Lowest NVR in a Nonwoven Blended Wiper

Product	Item	Description	Pack
55% Cell	TX609	TechniCloth™, 9x9	300 wipes/bag, 10 bags/cs
45% Poly	TX612	TechniCloth™, 12x12	150 wipes/bag, 10 bags/cs

DiNovo is an exclusive ITW Texwipe Distributor and sell an array of dry wipers. Contact us for all special requests.



www.agapilabsolutions.com

Contact Us:

Phone 1 (267) 397-3345

E-mail customerservice@agapilabsolutions.com

Wet Cleanroom Wipers



TechniSat®

TechniCloth® Wipers, Lowest NVR in a Nonwoven Blend, Pre-Wetted with Isopropyl Alcohol/DIW Solutions

Product	Item	Description	Pack
Non-Sterile	TX1065	9" x 11", Pre-wet with 70% IPA	50 wipes/flexpack; 24 flexpacks/cs
	TX1067	9" x 11", Pre-wet with 70% IPA	200 wipes/slider bg, 4 slider bg/cs

DiNovo is an exclusive ITW Texwipe Distributor and sell an array of wet wipers. Contact us for all special requests.



Sterile TechniSat®

TechniCloth® Wipers Pre-Wetted with 70% Isopropyl Alcohol and Sterilized

Product	Item	Description	Pack
Sterile	TX3214	9" x 11", Pre-wet with 70% IPA, STERILIZED	50 wipes/flexpack; 24 flexpacks/cs

DiNovo is an exclusive ITW Texwipe Distributor and sell an array of wet wipers. Contact us for all special requests.

Description

TechniSat® incorporates a pre-wetted fabric blend of nonwoven hydroentangled 55% cellulose and 45% polyester. This fabric combines the highly absorbent properties of a natural fiber with the cleanliness and strength of a synthetic fiber. It offers low particle counts, low extractable levels and solvent compatibility. TechniSat® wipers are pre-wetted with isopropyl alcohol/deionized water solutions filtered to 0.2 µm. The wipers are packaged in resealable flexpacks, flip-top tubs and canisters. They can be ordered as non-sterile or sterile.

Features

- Wipers pre-wetted with accurately diluted, measured volumes of high-purity IPA/deionized water solution
- Convenient, resealable dispensers
- Packaged in SolventSafe plastic that can be wiped down prior to introduction into the cleanroom
- Variety of configurations to a wide range of applications

Benefits

- Eliminates problems associated with purchasing, storing, filtering, diluting and mixing alcohol
- Eliminates maintenance of squirt bottles
- Offers portion and VOC* control with optimally wetted wipers
- Delivers premeasured amount of high-purity alcohol to the surface
- Provides strong, soft pre-wetted wiper for rapid wipe downs
- Offers low particle generation and low levels of solvent extractables



www.agapilabsolutions.com

Contact Us:

Phone 1 (267) 397-3345

E-mail customerservice@agapilabsolutions.com


Tamper Evident Sterile Seals

Steri-Tamp™ is safe, sterile, tamper evident seals for admixture ports on I.V. bags, bottles and vials. Protects against contamination and accidental double dosing. Alerts staff when a drug has been added. Foil outer layer provides visible evidence of a needle puncture. Removing the Steri-Tamp™ outer layer leaves a printed “OPENED” warning and cannot be reapplied.



28 mm Vial Seal

DiNovo Part Number	DiNovo Description	IVA® Seal Cross Reference
ST2103	13 mm Vial Seal (Red)	CP3003 Series
ST900901	20 mm Vial Seal (Silver)	CP3002 Series
ST1102	28 mm Vial Seal (Blue)	CP3001 Series
ST5103	IV Bag Port Seal (Red)	CP3011 & CP3013 Series
ST5112	IV Bag Port Seal (Blue)	CP3011 & CP3013 Series
ST4103	Syringe Seal (Red)	CP3020 & CP3030 Series
ST6101	Belly Button Bag Seal (Green)	CP3004 Series
ST7103	IV Bag Port Seal Chemo (Yellow)	CP3111Y



agapi
LAB SOLUTIONS
www.agapilabsolutions.com

Contact Us:
Phone 1 (267) 397-3345
E-mail customerservice@agapilabsolutions.com

Tamper Evident Syringe Caps



Tamper Evident Caps for IV syringes and oral dispensers are designed to increase doctor/nurse confidence while reducing cost. Tampering is immediately evident by a broken or missing sleeve. Caps are easily dispensed in supportive container. Caps are latex-free.

- Makes tampering immediately evident when the colored outer sleeve is broken or missing
- Increases the confidence level for Health Care Practitioners
- Allows filled dispensers with caps intact to be returned to the pharmacy for recycling



Innovating Tamper Evident Solutions to Enhance USP <797> Compliance

Product	Item	Description	Pack
Tamper Evident Caps for IV Syringes	D57-24	Bulk, Luer Lock, Sterilized, Red, Hard Pack	10/pk, 100pk/cs
	D57-25	Bulk, Luer Lock, Sterilized, Blue, Hard Pack	10/pk, 100pk/cs
	D57-26	Bulk, Luer Lock, Sterilized, White, Hard Pack	10/pk, 100pk/cs
	D57-30	Bulk, Luer Lock, Sterilized, Black, Hard Pack	10/pk, 100pk/cs
	D57-80	Individual, Luer Lock, Sterilized, Red, Blister Pk	100/box
Tamper Evident Caps for Oral Dispensers	D57-70	Bulk, Oral, BAXA, Red, Zip Bag	100/bag
	D54-74	Bulk, Oral, COMAR, Red, Zip Bag	100/bag



www.agapilabsolutions.com


Contact Us:

Phone 1 (267) 397-3345


E-mail customerservice@agapilabsolutions.com

Custom Lotion Pump Shippers

DiNovo has collaborated with compounding pharmacies and cold chain experts alike to develop a cost efficient, polyethylene shipper to protect and ship your lotion pumps with confidence. This cooler insert is designed specifically for lotion pumps making void fill, such as peanuts or bubble wrap, no longer necessary. The three pump slots are outfitted with die cut inserts to accommodate 50ml-150ml pumps. Packing and shipping your compounded lotion prescriptions have never been easier!



DiNovo's Standard Pack Out



Benefits

- Low profile design to reduce outbound freight cost
- 3 pump slot design with adjustable inserts to accommodate 50ml-150ml lotion pumps (GAKO, MEGA, and Topi-CLICK)
- Space to pack 8oz gel pack to keep product cold during shipments
- Eliminates the need for void fill, as the product is extremely secured and protected
- Vacuum Sealed Inserts to reduce inbound freight and warehouse space
- Need to alter the design to fit your pharmacy's needs? DiNovo has the ability to provide a fully custom product specifically designed for your pharmacy including box color, labels, and dimensions.

Product	Item	Description	Pack
Lotion Pump Shipper	LVP3SC	Lotion Pack, 8x8x4, 3-Multi-Size Slots, Case Pack	25/cs
	LVP3SC-G	Lotion Pack, 8x8x4, 3-Multi-Size Slots, 25 * 8oz Gel Packs, Case Pack	25/cs

**** Case Pack includes 25 individual coolers with all the components shown in diagram above. ****



www.agapilabsolutions.com


Contact Us:

Phone 1 (267) 397-3345

E-mail customerservice@agapilabsolutions.com

Packaging Products

- Corrugate
- EPS Products
- Polyethylene Products
- Gel Packs
- Molded Coolers
- Labels
- Bubble Wrap
- Tape
- Peanuts
- Plastic Film
- Zip Lock Bags



DiNovo, partnered with Houston Foam and Plastics, have the ability to fully manufacture any corrugate, EPS, and Polyethylene product. We do carry many of these items stock in our warehouse such as general size boxes and paneled coolers, but our ability to provide the perfect product for you is endless.

Product	Description
Corrugated Boxes	Full Manufacturing Ability: Size, Color, Print, Weight
EPS	Full Manufacturing Ability: Size, Thickness, Density, Design
Polyethylene	Full Manufacturing Ability: Size, Thickness, Density, Design



Lifoam Industries is the leader in custom foam molding and design engineering services for 50 years. Whether your need is protective packaging, temperature sensitive shipping, refrigerants, or custom foam components, DiNovo can distribute the product you need at the low cost you desire.

Product	Description
Gel Packs	Available Sizes: 8oz, 12oz, 24oz, 48oz
Molded Coolers	Array of .5"-3" Molded Coolers Delivered With or Without Box



DiNovo, partnered with Label Products, have a large array of stock labels that can be used for your pharmacy production or shipping needs. Need a custom label? No problem. DiNovo has the ability to produce your product in a week.

Product	Description
Labels	Full Design Process: Size, Color, Paper Stock, Finish, Text, Design



DiNovo, partnered with Ernest Packaging, has an extensive product catalog addressing virtually every packaging need. Bubble wrap, packaging tape, peanuts, stretch wrap, and so many more products are available in stock and ready to be shipped to you today.

Product	Description
Bubble Wrap	Available Sizes: 1/2" and 3/16" Bubble, Slit and Perforated
Tape	Available Tape: Hot Melt and Acrylic Tape, MIL Varies
Peanuts	Available Product: 20 cu.ft. Anti-Static, Green/Pink Peanuts
Stretch Wrap	Available Plastic Wrap: 18" x 1,500', 60GA, 70GA, 80GA, 90GA



DiNovo is a full service distributor of Elkay Plastics whom has the widest selection and broadest line of flexible plastic packaging in stock and available for immediate shipment. These products are designed for the healthcare marketplace meeting your needs for packaging. Need a custom bag? DiNovo can design and stock your specific bag needed.

Product	Description
Zip Lock Bags	Full Manufacturing Ability: Size, Density, Print



www.agapilabsolutions.com

Contact Us:

Phone 1 (267) 397-3345

E-mail customerservice@agapilabsolutions.com

PP DiNovo

Pharmacy & Packaging Provisions



www.agapilabsolutions.com

Contact Us:

Phone 1 (267) 397-3345

E-mail customerservice@agapilabsolutions.com