



# NEW! Locking Microcentrifuge Tubes

Available In 0.5mL, 1.5mL, 2.0mL & 5mL MACRO Tube

New Globe Scientific Locking Microcentrifuge Tubes are engineered with a secure locking system that provides the security of a screw cap microtube with the simplicity of a snap cap tube.



## Features at a Glance:

**Secure** – Maintains sample integrity with locking cap which prevents accidental cap opening and evaporation during centrifugation

**Versatile** – Single use tubes for laboratory procedures including centrifugation, boiling, incubation, freezing and sample storage

**Quality** – Tubes are made from optimally clear USP Class VI medical grade polypropylene

**Certified** – Manufactured in a class 100,000 clean room and are certified to be free from RNase and DNase, Human DNA, Pyrogens, PCR Inhibitors, and ATP

**Ease of Use** – Tubes feature graduations, writing area, pierceable cap, are autoclavable, and are packaged in zip, re-closable bags

Item #	Description	Unit
111552LK	0.5mL, polypropylene (PP), natural,	500
111562LK	1.5mL, polypropylene (PP), natural	500
111572LK	2.0mL, polypropylene (PP), natural	500
111582LK	5.0mL, polypropylene (PP), natural	100

### Contact Us:

Phone 1 (267) 397-3345

E-mail [customerservice@agopilabsolutions.com](mailto:customerservice@agopilabsolutions.com)